

April 9, 2010

To Whom It May Concern:

This is to hereby certify that the SeQual Technologies Inc. Eclipse Power Cartridge complies with IATA Guidance Document dated 10 August 2009 and titled "Transport of Lithium Metal and Lithium Ion Batteries.

The Eclipse Power Cartridge contains a quantity of two (2) independent batteries that are electrically separated until placed in the Eclipse system. Each battery contains 7.92 grams of equivalent Lithium Ion content. Each battery is 97.5 watt/hour which is below the guidance document limit of 100 watt/hour.

Brian Jarkell

Director, Quality Assurance & Regulatory Affairs

SeQual Technologies Inc